




Phospho-TNS1 (Y1326) Antibody

Product Code	CSB-PA000750
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HBL0
Immunogen	Synthesized peptide derived from Human Tensin-1 around the phosphorylation site of Y1326.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TNS1; TNS; Tensin-1
Immunogen Species	Homo sapiens (Human)
Target Names	TNS1

Image



Western Blot analysis of K562, cells using Phospho-Tensin-1 (Y1326) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.