




Phospho-ANXA1 (Y21) Antibody

Product Code	CSB-PA000753
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04083
Immunogen	Synthesized peptide derived from Human Annexin I around the phosphorylation site of Y21.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ANXA1; ANX1; LPC1; Annexin A1; Annexin I; Annexin-1; Calpactin II; Calpactin-2; Chromobindin-9; Lipocortin I; Phospholipase A2 inhibitory protein; p35
Immunogen Species	Homo sapiens (Human)
Target Names	ANXA1

Image



Western Blot analysis of NIH-3T3, cells using Phospho-Annexin I (Y21) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.