




Phospho-MAP3K11 (S674) Antibody

Product Code	CSB-PA000754
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16584
Immunogen	Synthesized peptide derived from Human MLK3 around the phosphorylation site of S674.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MAP3K11; MLK3; PTK1; SPRK; Mitogen-activated protein kinase kinase kinase 11; Mixed lineage kinase 3; Src-homology 3 domain-containing proline-rich kinase
Immunogen Species	Homo sapiens (Human)
Target Names	MAP3K11

Image



Western Blot analysis of NIH-3T3, cells using Phospho-MLK3 (S674) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.