




Phospho-CTNNB1 (S552) Antibody

Product Code	CSB-PA000756
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35222
Immunogen	Synthesized peptide derived from Human Catenin-β around the phosphorylation site of S552.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CTNNB1; CTNNB; OK/SW-cl.35; Catenin beta-1; Beta-catenin
Immunogen Species	Homo sapiens (Human)
Target Names	CTNNB1

Image



Western Blot analysis of NIH-3T3, cells using Phospho-Catenin-β (S552) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.