




Phospho-EPOR (Y426) Antibody

Product Code	CSB-PA000757
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19235
Immunogen	Synthesized peptide derived from Human EpoR around the phosphorylation site of Y426.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EPOR; Erythropoietin receptor; EPO-R
Immunogen Species	Homo sapiens (Human)
Target Names	EPOR

Image



Western Blot analysis of K562, cells using Phospho-EpoR (Y426) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.