





HSD3B7 Antibody

Product Code	CSB-PA000767
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2F3
Immunogen	Synthesized peptide derived from the Internal region of Human 3 β -HSD7.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HSD3B7

Image



Western Blot analysis of HeLa cells using 3 β -HSD7 Polyclonal Antibody



Western Blot analysis of MCF7 HELA cells using 3 β -HSD7 Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.