





EPB41L2 Antibody

Product Code	CSB-PA000768
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43491
Immunogen	Synthesized peptide derived from the Internal region of Human 4.1G.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EPB41L2; Band 4.1-like protein 2; Generally expressed protein 4.1; 4.1G
Immunogen Species	Homo sapiens (Human)
Target Names	EPB41L2

Image



Western Blot analysis of K562 cells using 4.1G Polyclonal Antibody



Western Blot analysis of K562 cells using 4.1G Polyclonal Antibody

