





EPB41L2 Antibody

Product Code	CSB-PA000768
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43491
Immunogen	Synthesized peptide derived from the Internal region of Human 4.1G.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	EPB41L2; Band 4.1-like protein 2; Generally expressed protein 4.1; 4.1G
Immunogen Species	Homo sapiens (Human)
Target Names	EPB41L2

Image



Western Blot analysis of K562 cells using 4.1G Polyclonal Antibody



Western Blot analysis of K562 cells using 4.1G Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.