




TP53BP1 Antibody

Product Code	CSB-PA000771
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12888
Immunogen	Synthesized peptide derived from Human 53BP1 around the non-phosphorylation site of S6.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TP53BP1; Tumor suppressor p53-binding protein 1; 53BP1; p53-binding protein 1; p53BP1
Immunogen Species	Homo sapiens (Human)
Target Names	TP53BP1

Image



Western Blot analysis of A549 cells using 53BP1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.