





TNFAIP3 Antibody

Product Code	CSB-PA000775
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21580
Immunogen	Synthesized peptide derived from the Internal region of Human A20.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TNFAIP3; OTUD7C; Tumor necrosis factor alpha-induced protein 3; TNF alpha-induced protein 3; OTU domain-containing protein 7C; Putative DNA-binding protein A20; Zinc finger protein A20
Immunogen Species	Homo sapiens (Human)
Target Names	TNFAIP3

Image



Western Blot analysis of 293 cells using A20 Polyclonal Antibody



Western Blot analysis of various cells using A20 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.