




AAK1 Antibody

Product Code	CSB-PA000776
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2M2I8
Immunogen	Synthesized peptide derived from the Internal region of Human AAK1.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Monkey
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	AAK1; KIAA1048; AP2-associated protein kinase 1; Adaptor-associated kinase 1
Immunogen Species	Homo sapiens (Human)
Target Names	AAK1

Image



Western Blot analysis of COS7 cells using AAK1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.