




ABCA6 Antibody

Product Code	CSB-PA000780
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N139
Immunogen	Synthesized peptide derived from the C-terminal region of Human ABCA6.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ABCA6; ATP-binding cassette sub-family A member 6
Immunogen Species	Homo sapiens (Human)
Target Names	ABCA6

Image



Western Blot analysis of 293 cells using ABCA6 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.