





ABCA8 Antibody

Product Code	CSB-PA000781
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94911
Immunogen	Synthesized peptide derived from the Internal region of Human ABCA8.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ABCA8; KIAA0822; ATP-binding cassette sub-family A member 8
Immunogen Species	Homo sapiens (Human)
Target Names	ABCA8

Image



Western Blot analysis of HuvEc cells using ABCA8 Polyclonal Antibody



Western Blot analysis of HuvEc PC12 cells using ABCA8 Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.