





ABCC13 Antibody

Product Code	CSB-PA000785
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NSE7
Immunogen	Synthesized peptide derived from the Internal region of Human ABCC13.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ABCC13; C21orf73; PRED6; Putative ATP-binding cassette sub-family C member 13
Immunogen Species	Homo sapiens (Human)
Target Names	ABCC13

Image



Western Blot analysis of HT29 cells using ABCC13 Polyclonal Antibody



Western Blot analysis of HELA cells using ABCC13 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.