





ABCD4 Antibody

Product Code	CSB-PA000786
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14678
Immunogen	Synthesized peptide derived from the Internal region of Human ABCD4.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ABCD4; PXMP1L; ATP-binding cassette sub-family D member 4; PMP70-related protein; P70R; Peroxisomal membrane protein 1-like; PXMP1-L; Peroxisomal membrane protein 69; PMP69
Immunogen Species	Homo sapiens (Human)
Target Names	ABCD4

Image



Western Blot analysis of HepG2 cells using ABCD4 Polyclonal Antibody



Western Blot analysis of K562 cells using
ABCD4 Polyclonal Antibody