




ABHD11 Antibody

Product Code	CSB-PA000789
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFV4
Immunogen	Synthesized peptide derived from the Internal region of Human ABHD11.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ABHD11; WBSCR21; PP1226; Abhydrolase domain-containing protein 11; Williams-Beuren syndrome chromosomal region 21 protein
Immunogen Species	Homo sapiens (Human)
Target Names	ABHD11

Image



Western Blot analysis of HuvEc cells using ABHD11 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.