




ABHD12 Antibody

Product Code	CSB-PA000790
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N2K0
Immunogen	Synthesized peptide derived from the Internal region of Human ABHD12.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ABHD12; C20orf22; Monoacylglycerol lipase ABHD12; 2-arachidonoylglycerol hydrolase; Abhydrolase domain-containing protein 12
Immunogen Species	Homo sapiens (Human)
Target Names	ABHD12

Image



Western Blot analysis of COLO205 cells using ABHD12 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.