




DBI Antibody

Product Code	CSB-PA000798
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07108
Immunogen	Synthesized peptide derived from the Internal region of Human ACBP.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DBI; Acyl-CoA-binding protein; ACBP; Diazepam-binding inhibitor; DBI; Endozepine; EP
Immunogen Species	Homo sapiens (Human)
Target Names	DBI

Image



Western Blot analysis of HuvEc cells using ACBP Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.