





CHRNA10 Antibody

Product Code	CSB-PA000800
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZZ6
Immunogen	Synthesized peptide derived from the C-terminal region of Human AChR α 10.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CHRNA10; NACHRA10; Neuronal acetylcholine receptor subunit alpha-10; Nicotinic acetylcholine receptor subunit alpha-10; NACHR alpha-10
Immunogen Species	Homo sapiens (Human)
Target Names	CHRNA10

Image



Western Blot analysis of A549 cells using AChR α 10 Polyclonal Antibody



Western Blot analysis of HuvEc MCF7 cells using AChR α 10 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.