





ADIPOQ Antibody

Product Code	CSB-PA000806
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15848
Immunogen	Synthesized peptide derived from the Internal region of Human Acrp30.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ADIPOQ; ACDC; ACRP30; APM1; GBP28; Adiponectin; 30 kDa adipocyte complement-related protein; Adipocyte complement-related 30 kDa protein; ACRP30; Adipocyte; C1q and collagen domain-containing protein; Adipose most abundant gene transcript 1
Immunogen Species	Homo sapiens (Human)
Target Names	ADIPOQ

Image



Western Blot analysis of HeLa cells using Acrp30 Polyclonal Antibody



Western Blot analysis of VEC 293 cells using Acrp30 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.