





ACTA1 Antibody

Product Code	CSB-PA000810
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68133
Immunogen	Synthesized peptide derived from the N-terminal region of Human Actin α 1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ACTA1; ACTA; Actin; alpha skeletal muscle; Alpha-actin-1
Immunogen Species	Homo sapiens (Human)
Target Names	ACTA1

Image



Western Blot analysis of HEPG2-UV cells using Actin α 1 Polyclonal Antibody



Western Blot analysis of LOVO mouse-muscle cells using Actin α 1 Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.