





ACTB Antibody

Product Code	CSB-PA000813
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60709
Immunogen	Synthesized peptide derived from the N-terminal region of Human Actin β .
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:1000-1:4000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ACTB; Actin; cytoplasmic 1; Beta-actin
Immunogen Species	Homo sapiens (Human)
Target Names	ACTB


Image



Western Blot analysis of HeLa cells using Actin β Polyclonal Antibody



Western Blot analysis of HELA cells using Actin β Polyclonal Antibody



Western blot analysis of customer