





ACTN1/ACTN2/ACTN3/ACTN4 Antibody

Product Code	CSB-PA000814
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12814/P35609/Q08043/O43707
Immunogen	Synthesized peptide derived from the N-terminal region of Human Actinin- α 1/2/3/4.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ACTN1; Alpha-actinin-1; Alpha-actinin cytoskeletal isoform; F-actin cross-linking protein; Non-muscle alpha-actinin-1; ACTN2; Alpha-actinin-2; Alpha-actinin skeletal muscle isoform 2; F-actin cross-linking protein; ACTN3; Alpha-actinin-3; A
Immunogen Species	Homo sapiens (Human)
Target Names	ACTN1/ACTN2/ACTN3/ACTN4

Image



Western Blot analysis of 293 cells using Actinin- α 1/2/3/4 Polyclonal Antibody



Western Blot analysis of Jurkat 293 HELA A549 cells using Actinin- α 1/2/3/4 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.