




ADAM32 Antibody

Product Code	CSB-PA000819
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TC27
Immunogen	Synthesized peptide derived from the Internal region of Human ADAM32.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ADAM32; Disintegrin and metalloproteinase domain-containing protein 32; ADAM 32
Immunogen Species	Homo sapiens (Human)
Target Names	ADAM32

Image



Western Blot analysis of 293 cells using ADAM32 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.