





AHR/AHRR Antibody

Product Code	CSB-PA000832
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35869/A9YTQ3
Immunogen	Synthesized peptide derived from Human Ah Receptor around the non-phosphorylation site of S36.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	AHR; BHLHE76; Aryl hydrocarbon receptor; Ah receptor; AhR; Class E basic helix-loop-helix protein 76; bHLHe76; AHRR; BHLHE77; KIAA1234; Aryl hydrocarbon receptor repressor; AhR repressor; AhRR; Class E basic helix-loop-helix protein 77; bHL
Immunogen Species	Homo sapiens (Human)
Target Names	AHR/AHRR

Image



Western Blot analysis of 293 cells using Ah Receptor Polyclonal Antibody



Western Blot analysis of 293 cells using Ah Receptor Polyclonal Antibody