




AK5 Antibody

Product Code	CSB-PA000840
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6K8
Immunogen	Synthesized peptide derived from the C-terminal region of Human AK5.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	AK5; Adenylate kinase isoenzyme 5; AK 5; ATP-AMP transphosphorylase 5
Immunogen Species	Homo sapiens (Human)
Target Names	AK5

Image



Western Blot analysis of Jurkat cells using AK5 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.