




AKAP5 Antibody

Product Code	CSB-PA000843
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24588
Immunogen	Synthesized peptide derived from the N-terminal region of Human AKAP 79.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	AKAP5; AKAP79; A-kinase anchor protein 5; AKAP-5; A-kinase anchor protein 79 kDa; AKAP 79; H21; cAMP-dependent protein kinase regulatory subunit II high affinity-binding protein
Immunogen Species	Homo sapiens (Human)
Target Names	AKAP5

Image



Western Blot analysis of COLO205 cells using AKAP 79 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.