




AKR1E2 Antibody

Product Code	CSB-PA000847
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96JD6
Immunogen	Synthesized peptide derived from the Internal region of Human AKR1CL2.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	AKR1E2; AKR1CL2; AKRDC1; 1; 5-anhydro-D-fructose reductase; AF reductase; Aldo-keto reductase family 1 member C-like protein 2; Aldo-keto reductase family 1 member E2; LoopADR; Testis-specific protein; hTSP
Immunogen Species	Homo sapiens (Human)
Target Names	AKR1E2

Image



Western Blot analysis of A549 cells using AKR1CL2 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.