





AKT1 Antibody

Product Code	CSB-PA000849
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31749
Immunogen	Synthesized peptide derived from Human Akt1 around the non-phosphorylation site of S246.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AKT1


Image




Western Blot analysis of 293 cells using Akt1 Polyclonal Antibody




Western Blot analysis of Jurkat cells using Akt1 Polyclonal Antibody



Western Blot analysis of HELA KB PC-3 using Akt1 Polyclonal Antibody.



Western Blot analysis of HELA HEPG2 PC-3 using Akt1 Polyclonal Antibody.



Western Blot analysis of HEPG2 using Akt1 Polyclonal Antibody