





AKT1 Antibody

Product Code	CSB-PA000851
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31749
Immunogen	Synthesized peptide derived from Human Akt around the non-phosphorylation site of Y326.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Monkey
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	AKT1; PKB; RAC; RAC-alpha serine/threonine-protein kinase; Protein kinase B; PKB; Protein kinase B alpha; PKB alpha; Proto-oncogene c-Akt; RAC-PK-alpha
Immunogen Species	Homo sapiens (Human)
Target Names	AKT1

Image



Western Blot analysis of COS7 cells using Akt Polyclonal Antibody



Western Blot analysis of HEPG2 using Akt Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.