





ALDH1B1 Antibody

Product Code	CSB-PA000860
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30837
Immunogen	Synthesized peptide derived from the Internal region of Human ALDH1B1.
Raised In	Rabbit
Species Reactivity	Human,Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ALDH1B1; ALDH5; ALDHX; Aldehyde dehydrogenase X; mitochondrial; Aldehyde dehydrogenase 5; Aldehyde dehydrogenase family 1 member B1
Immunogen Species	Homo sapiens (Human)
Target Names	ALDH1B1

Image



Western Blot analysis of COS7 cells using ALDH1B1 Polyclonal Antibody



Western Blot analysis of HELA cells using ALDH1B1 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.