




ANKRD30A Antibody

Product Code	CSB-PA000876
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXX3
Immunogen	Synthesized peptide derived from the C-terminal region of Human ANKRD30A.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ANKRD30A; Ankyrin repeat domain-containing protein 30A; Serologically defined breast cancer antigen NY-BR-1
Immunogen Species	Homo sapiens (Human)
Target Names	ANKRD30A

Image



Western Blot analysis of COLO205 cells using ANKRD30A Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.