




NUSAP1 Antibody

Product Code	CSB-PA000877
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXS6
Immunogen	Synthesized peptide derived from the C-terminal region of Human ANKT.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NUSAP1; ANKT; BM-037; Nucleolar and spindle-associated protein 1; NuSAP
Immunogen Species	Homo sapiens (Human)
Target Names	NUSAP1

Image



Western Blot analysis of HT29 cells using ANKT Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.