





ANXA6 Antibody

Product Code	CSB-PA000881
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08133
Immunogen	Synthesized peptide derived from the N-terminal region of Human Annexin VI.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ANXA6; ANX6; Annexin A6; 67 kDa calelectrin; Annexin VI; Annexin-6; Calphobindin-II; CPB-II; Chromobindin-20; Lipocortin VI; Protein III; p68; p70
Immunogen Species	Homo sapiens (Human)
Target Names	ANXA6

Image



Western Blot analysis of 293 cells using Annexin VI Polyclonal Antibody



Western Blot analysis of 293 cells using Annexin VI Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.