




TFAP2C Antibody

Product Code	CSB-PA000894
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92754
Immunogen	Synthesized peptide derived from the C-terminal region of Human AP-2 γ .
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TFAP2C; Transcription factor AP-2 gamma; AP2-gamma; Activating enhancer-binding protein 2 gamma; Transcription factor ERF-1
Immunogen Species	Homo sapiens (Human)
Target Names	TFAP2C

Image



Western Blot analysis of HuvEc cells using AP-2 γ Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.