





APAF1 Antibody

Product Code	CSB-PA000895
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14727-6
Immunogen	Synthesized peptide derived from the C-terminal region of Human Apaf-1-ALT.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	APAF1-interacting protein; APIP; APIP2; CGI-29 protein; MMRP19; dJ179L10.2; likely ortholog of mouse monocyte macrophage 19
Immunogen Species	Homo sapiens (Human)
Target Names	APAF1

Image



Western Blot analysis of COLO205 cells using Apaf-1-ALT Polyclonal Antibody



Western Blot analysis of Jurkat cells using Apaf-1-ALT Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.