





APC Antibody

Product Code	CSB-PA000897
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25054
Immunogen	Synthesized peptide derived from the C-terminal region of Human APC.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	APC; DP2.5; Adenomatous polyposis coli protein; Protein APC; Deleted in polyposis 2.5
Immunogen Species	Homo sapiens (Human)
Target Names	APC

Image



Western Blot analysis of COLO205 cells using APC Polyclonal Antibody



Western Blot analysis of RAT-MUSCLE cells using APC Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.