





APOF Antibody

Product Code	CSB-PA000905
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13790
Immunogen	Synthesized peptide derived from the C-terminal region of Human ApoF.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	APOF; Apolipoprotein F; Apo-F; Lipid transfer inhibitor protein; LTIP
Immunogen Species	Homo sapiens (Human)
Target Names	APOF

Image



Western Blot analysis of A549 cells using ApoF Polyclonal Antibody



Western Blot analysis of KB PC-3 MCF7 cells using ApoF Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.