




APPL1 Antibody

Product Code	CSB-PA000908
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKG1
Immunogen	Synthesized peptide derived from the Internal region of Human APPL1.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	APPL1; APPL; DIP13A; KIAA1428; DCC-interacting protein 13-alpha; Dip13-alpha; Adapter protein containing PH domain; PTB domain and leucine zipper motif 1
Immunogen Species	Homo sapiens (Human)
Target Names	APPL1

Image



Western Blot analysis of A549 cells using APPL1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.