





TNFSF13 Antibody

Product Code	CSB-PA000909
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75888
Immunogen	Synthesized peptide derived from the Internal region of Human APRIL.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TNFSF13; APRIL; TALL2; ZTNF2; Tumor necrosis factor ligand superfamily member 13; A proliferation-inducing ligand; APRIL; TNF- and APOL-related leukocyte expressed ligand 2; TALL-2; TNF-related death ligand 1; TRDL-1; CD antigen CD256
Immunogen Species	Homo sapiens (Human)
Target Names	TNFSF13

Image



Western Blot analysis of COLO205 cells using APRIL Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using APRIL Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.