




MIP Antibody

Product Code	CSB-PA000910
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30301
Immunogen	Synthesized peptide derived from the Internal region of Human AQP0.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MIP; AQP0; Lens fiber major intrinsic protein; Aquaporin-0; MIP26; MP26
Immunogen Species	Homo sapiens (Human)
Target Names	MIP

Image



Western Blot analysis of HT29 cells using AQP0 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.