




AQP12A/AQP12B Antibody

Product Code	CSB-PA000912
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IXF9/A6NM10
Immunogen	Synthesized peptide derived from the C-terminal region of Human AQP12.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	AQP12A; AQP12; AQPX2; Aquaporin-12A; AQP-12; AQP12B; Aquaporin-12B; AQP-12B
Immunogen Species	Homo sapiens (Human)
Target Names	AQP12A/AQP12B

Image



Western Blot analysis of Jurkat cells using AQP12 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.