





ARHGAP18 Antibody

Product Code	CSB-PA000921
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N392
Immunogen	Synthesized peptide derived from the C-terminal region of Human ARHGAP18.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ARHGAP18; Rho GTPase-activating protein 18; MacGAP; Rho-type GTPase-activating protein 18
Immunogen Species	Homo sapiens (Human)
Target Names	ARHGAP18

Image



Western Blot analysis of Jurkat cells using ARHGAP18 Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using ARHGAP18 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.