




ADRB2 Antibody

Product Code	CSB-PA000938
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07550
Immunogen	Synthesized peptide derived from Human AR-β2 around the non-phosphorylation site of S346.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, FCM, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, FC:1:200-1:400, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ADRB2; ADRB2R; B2AR; Beta-2 adrenergic receptor; Beta-2 adrenoreceptor; Beta-2 adrenoceptor
Immunogen Species	Homo sapiens (Human)
Target Names	ADRB2

Image



Western Blot analysis of BxPc3 cells using AR-β2 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.