





MAP3K5 Antibody

Product Code	CSB-PA000942
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99683
Immunogen	Synthesized peptide derived from Human ASK 1 around the non-phosphorylation site of S83.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MAP3K5; ASK1; MAPKKK5; MEKK5; Mitogen-activated protein kinase kinase 5; Apoptosis signal-regulating kinase 1; ASK-1; MAPK/ERK kinase 5; MEK kinase 5; MEKK 5
Immunogen Species	Homo sapiens (Human)
Target Names	MAP3K5

Image



Western Blot analysis of 293T cells using ASK 1 Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using ASK 1 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.