




ATPAF2 Antibody

Product Code	CSB-PA000956
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N5M1
Immunogen	Synthesized peptide derived from the Internal region of Human ATPAF2.
Raised In	Rabbit
Species Reactivity	Human,Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ATPAF2; ATP12; LP3663; ATP synthase mitochondrial F1 complex assembly factor 2; ATP12 homolog
Immunogen Species	Homo sapiens (Human)
Target Names	ATPAF2

Image



Western Blot analysis of Jurkat cells using ATPAF2 Polyclonal Antibody