




AVEN Antibody

Product Code	CSB-PA000961
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQS1
Immunogen	Synthesized peptide derived from the C-terminal region of Human Aven.
Raised In	Rabbit
Species Reactivity	Human,Monkey
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	AVEN; Cell death regulator Aven
Immunogen Species	Homo sapiens (Human)
Target Names	AVEN

Image



Western Blot analysis of LOVO cells using Aven Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.