




BACH1 Antibody

Product Code	CSB-PA000965
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14867
Immunogen	Synthesized peptide derived from the Internal region of Human BACH1.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BACH1; Transcription regulator protein BACH1; BTB and CNC homolog 1; HA2303
Immunogen Species	Homo sapiens (Human)
Target Names	BACH1

Image



Western Blot analysis of 293 cells using BACH1 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.