





# BAK1 Antibody

<b>Product Code</b>	CSB-PA000971
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q16611
<b>Immunogen</b>	Synthesized peptide derived from the N-terminal region of Human Bak.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	BAK1; BAK; BCL2L7; CDN1; Bcl-2 homologous antagonist/killer; Apoptosis regulator BAK; Bcl-2-like protein 7; Bcl2-L-7
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	BAK1

## Image



Western Blot analysis of 293 cells using Bak Polyclonal Antibody



Western Blot analysis of 293T HELA MOUSE-BRAIN cells using Bak Polyclonal Antibody