





BAK1 Antibody

Product Code	CSB-PA000971
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16611
Immunogen	Synthesized peptide derived from the N-terminal region of Human Bak.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BAK1; BAK; BCL2L7; CDN1; Bcl-2 homologous antagonist/killer; Apoptosis regulator BAK; Bcl-2-like protein 7; Bcl2-L-7
Immunogen Species	Homo sapiens (Human)
Target Names	BAK1

Image



Western Blot analysis of 293 cells using Bak Polyclonal Antibody



Western Blot analysis of 293T HELA MOUSE-BRAIN cells using Bak Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.